

Anti-CD3 antibody [17A2] {PE/Cyanine7} (STJA0000186)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine7
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12503 12501 12500 Cd247 Cd3e CD3Z_MOUSE CD3E_MOUSE CD3D_MOUSE
---------	--

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---



C57BL/6 murine splenocytes were stained with Anti-CD3 antibody (STJA0000186) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.