

Anti-CD3 antibody [17A2] {PerCP/Cyanine5.5} (STJA0000187)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID
[12502](#)
[12500](#)
[12501](#)
[Cd3g](#)
[Cd3d](#)
[CD3G_MOUSE](#)
[CD3D_MOUSE](#)
[CD3E_MOUSE](#)

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine splenocytes were stained with Anti-CD3 antibody (STJA0000187) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.