

Anti-CD3 antibody [17A2] {AF647} (STJA0000198)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	AF647
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID
[12502](#)
[12500](#)
[12503](#)
[Cd3g](#)
[Cd3d](#)
[CD3G_MOUSE](#)
[CD3D_MOUSE](#)
[CD3Z_MOUSE](#)

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine splenocytes were stained with Anti-Mouse CD3 Monoclonal Antibody (AF647 Conjugated)[Used at 0.2 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained lymphocytes (empty black histogram) were used as control.