

## Anti-CD47 antibody [MIAP410] {PerCP} (STJA0000202)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PerCP
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">16423</a> <a href="#">29364</a> <a href="#">961</a> <a href="#">Cd47</a> <a href="#">Cd47</a> <a href="#">CD47_MOUSE</a> <a href="#">CD47_RAT</a> <a href="#">CD47_HUMAN</a>
---------	---

### ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---



C57BL/6 murine splenocytes were stained with Anti-CD47 antibody (STJA0000202) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.