

Anti-CD47 antibody [MIAP410] {AF488} (STJA0000213)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	AF488
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	961 16423 29364 CD47 Cd47 CD47_HUMAN CD47_MOUSE CD47_RAT
---------	---

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---



C57BL/6 murine splenocytes were stained with Anti-CD47 antibody (STJA0000213) [Used at 0.2 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.