

Anti-CD54/Icam1 antibody [YN1/1.7.4] {AF647} (STJA0000235)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	AF647
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	15894
Gene Symbol	Icam1
UniProt ID	ICAM1_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine splenocytes were stained with Anti-Mouse CD54 Monoclonal Antibody (AF647 Conjugated) (filled gray histogram). Unstained lymphocytes (empty black histogram) were used as control.