

Anti-CD127/IL7r antibody [A7R34] {AF647} (STJA0000260)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	AF647
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16197
Gene Symbol	IL7r
UniProt ID	IL7RA_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine splenocytes were stained with Anti-Mouse CD127 Monoclonal Antibody (AF647 Conjugated) [Used at 0.2 µg/10⁶ cells dilution] and Anti-Mouse CD3 Monoclonal Antibody (AF488 Conjugated). Splenocytes stained with Anti-Mouse CD3 Monoclonal Antibody (AF488 Conjugated) (right) were used as control.