

Anti-CD21/CR2 antibody [BU32] {FITC} (STJA0000332)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	1380
Gene Symbol	CR2
UniProt ID	CR2_HUMAN

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Human peripheral blood lymphocytes were stained with Anti-CD21 antibody (STJA0000332) (filled gray histogram). Unstained lymphocytes (empty black histogram) were used as control.