

Anti-CD15/FUT4 antibody [HI98] {APC} (STJA0000411)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	APC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Human peripheral blood neutrophils were stained with Anti-CD15 antibody (STJA0000411) (filled gray histogram). Unstained neutrophils (empty black histogram) were used as control.