

Anti-CD11b antibody [M1/70] {PE} (STJA0000422)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse/Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3684 NA
Gene Symbol	ITGAM Itgam
UniProt ID	ITAM_HUMAN ITAM_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Mouse bone marrow cells were stained with Anti-Mouse/Human CD11b antibody (STJA0000422) (filled gray histogram). Unstained bone marrow cells (blank black histogram) were used as control.