

## Anti-CD11b antibody [M1/70] {PerCP} (STJA0000424)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse/Human

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PerCP
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene Symbol	<a href="#">Itgam</a> <a href="#">ITGAM</a>
UniProt ID	<a href="#">ITAM_MOUSE</a> <a href="#">ITAM_HUMAN</a>

### ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Mouse bone marrow cells were stained with Anti-Mouse/Human CD11b Monoclonal Antibody (PerCP Conjugated) (filled gray histogram). Unstained bone marrow cells (blank black histogram) were used as control.