

Anti-CD11b antibody [M1/70] {PE/Cyanine5} (STJA0000434)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3684 NA
Gene Symbol	ITGAM Itgam
UniProt ID	ITAM_HUMAN ITAM_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Mouse bone marrow cells were stained with Anti-Mouse/Human CD11b antibody (STJA0000434) [Used at 0.01 $\mu\text{g}/10^6$ cells dilution] (filled gray histogram). Unstained bone marrow cells (blank black histogram) were used as control.