

## Anti-CD11b antibody [M1/70] {PE/Cyanine7} (STJA0000435)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Human/Mouse

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine7
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">3684</a> <a href="#">NA</a>
Gene Symbol	<a href="#">ITGAM</a> <a href="#">Itgam</a>
UniProt ID	<a href="#">ITAM_HUMAN</a> <a href="#">ITAM_MOUSE</a>

### ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine bone marrow cells were stained with Anti-Mouse/Human CD11b antibody (STJA0000435) [Used at 0.2  $\mu$ g/ $10^6$  cells dilution] (filled gray histogram). Unstained bone marrow cells (empty black histogram) were used as control.