

Anti-CD11b antibody [M1/70] {PerCP/Cyanine5.5} (STJA0000436)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PerCP/Cyanine5.5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene Symbol	Itgam ITGAM
UniProt ID	ITAM_MOUSE ITAM_HUMAN

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine bone marrow cells were stained with Anti-Mouse/Human CD11b antibody (STJA0000436) [Used at 0.05 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained bone marrow cells (empty black histogram) were used as control.