

## Anti-CD11b antibody [M1/70] {AF647} (STJA0000438)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Human/Mouse

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	AF647
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene Symbol	<a href="#">Itgam</a> <a href="#">ITGAM</a>
UniProt ID	<a href="#">ITAM_MOUSE</a> <a href="#">ITAM_HUMAN</a>

### ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine bone marrow cells were stained with Anti-Mouse/Human CD11b Monoclonal Antibody (AF647 Conjugated) [Used at 0.05  $\mu$ g/ $10^6$  cells dilution] (filled gray histogram). Unstained bone marrow cells (empty black histogram) were used as control.