

Anti-CD64 antibody [10.1] {PerCP/Cyanine5.5} (STJA0000443)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	2209 NA
Gene Symbol	FCGR1A FCGR1BP
UniProt ID	FCGR1_HUMAN FCGRB_HUMAN

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---



Human peripheral blood monocytes were stained with Anti-CD64 antibody (STJA0000443) (filled gray histogram). Unstained monocytes (empty black histogram) were used as control.