

Anti-CD103/Itgae antibody [M290] {AF647} (STJA0000462)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	AF647
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16407
Gene Symbol	Itgae
UniProt ID	ITAE_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine splenocytes were stained with Anti-Mouse CD103 Monoclonal Antibody (AF647 Conjugated) and Anti-Mouse CD3 Monoclonal Antibody (AF488 Conjugated) (left). Splenocytes stained with Anti-Mouse CD3 Monoclonal Antibody (AF488 Conjugated) (right) were used as control.