

Anti-CD117/Kit antibody [2B8] {PE} (STJA0000483)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16590
Gene Symbol	Kit
UniProt ID	KIT_MOUSE

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



C57BL/6 murine bone marrow cells were stained with Anti-CD117 antibody (STJA0000483) [Used at 0.2 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained bone marrow cells (empty black histogram) were used as control.