

## Anti-Cd4 antibody [GK1.5] {FITC} (STJA0000520)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rat                |
| Reactivity    | Mouse              |

### PRODUCT PROPERTIES


|                     |                                                                                                  |
|---------------------|--------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                       |
| Concentration       | 5 Mu L/Test                                                                                      |
| Conjugation         | FITC                                                                                             |
| Formulation         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.      |
| Isotype             | IgG2bk                                                                                           |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

### TARGET INFORMATION

|             |                           |
|-------------|---------------------------|
| Gene ID     | <a href="#">12504</a>     |
| Gene Symbol | <a href="#">Cd4</a>       |
| UniProt ID  | <a href="#">CD4_MOUSE</a> |

### ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Mouse splenocytes were stained with Anti-CD4 antibody (STJA0000520) (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.