

Anti-Cd4 antibody [GK1.5] {PE/Cyanine5} (STJA0000533)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12504
Gene Symbol	Cd4
UniProt ID	CD4_MOUSE

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Mouse splenocytes were stained with Anti-CD4 antibody (STJA0000533) [Used at 0.04 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.