

Anti-Cd4 antibody [GK1.5] {AF488} (STJA0000535)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	AF488
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12504
Gene Symbol	Cd4
UniProt ID	CD4_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Mouse splenocytes were stained with Anti-CD4 antibody (STJA0000535) [Used at 0.2 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.