

Anti-CD25/Il2ra antibody [PC-61.5.3] {FITC} (STJA0000568)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16184
Gene Symbol	Il2ra
UniProt ID	IL2RA_MOUSE

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



C57BL/6 murine splenocytes were stained with Anti-CD25 antibody (STJA0000568) [Used at 0.2 μ g/ 10^6 cells dilution] and Anti-Mouse CD4 Monoclonal Antibody (APC Conjugated) (right). Splenocytes stained with Anti-Mouse CD4 Monoclonal Antibody (APC Conjugated), are used as control.