

Anti-CD25/Il2ra antibody [PC-61.5.3] {APC} (STJA0000570)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	APC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16184
Gene Symbol	Il2ra
UniProt ID	IL2RA_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine splenocytes were stained with Anti-CD25 antibody (STJA0000570) [Used at 0.2 μ g/ 10^6 cells dilution] and Anti-Mouse CD4 Monoclonal Antibody (FITC Conjugated). Splenocytes stained with Anti-Mouse CD4 Monoclonal Antibody (FITC Conjugated) and Rat IgG1 Isotype Control (APC Conjugated) (right) were used as control.