

Anti-Cd8a antibody [53-6.7] {APC} (STJA0000598)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	APC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12525
Gene Symbol	Cd8a
UniProt ID	CD8A_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Mouse splenocytes were stained with Anti-CD8a antibody (STJA0000598) [Used at 0.01 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained splenocytes (blank black histogram) were used as control.