

Anti-Cd4 antibody [OX-38] {FITC} (STJA0000610)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	24932
Gene Symbol	Cd4
UniProt ID	CD4_RAT

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Rat splenocytes were stained with Anti-Rat CD4 (domain 1 Conjugated) Monoclonal Antibody (FITC Conjugated)[Used at 0.1 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.