

Anti-Ly6g antibody [1A8] {PE} (STJA0000624)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	546644
Gene Symbol	Ly6g
UniProt ID	LY6G_MOUSE

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---



Mouse bone marrow cells were stained with Anti-Ly6G antibody (STJA0000624) [Used at 0.02 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained bone marrow cells (blank black histogram) were used as control.