

Anti-Ly6g antibody [1A8] {PerCP} (STJA0000626)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PerCP
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	546644
Gene Symbol	Ly6g
UniProt ID	LY6G_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Mouse bone marrow cells were stained with Anti-Mouse Ly6G Monoclonal Antibody (PerCP Conjugated)[Used at 0.05 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained bone marrow cells (black histogram) were used as control.