

Anti-CD4 antibody [RPA-T4] {PE/Cyanine5.5} (STJA0000638)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES


| | |
|---------------------|--------------------------------------------------------------------------------------------------|
| Clonality | Monoclonal |
| Concentration | 5 Mu L/Test |
| Conjugation | PE/Cyanine5.5 |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG1k |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------|---------------------------|
| Gene ID | 920 |
| Gene Symbol | CD4 |
| UniProt ID | CD4_HUMAN |

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Human peripheral blood cells were stained with Anti-CD4 antibody (STJA0000638) (filled gray histogram). Unstained peripheral blood cells (blank black histogram) were used as control.