

Anti-CD45R/B220 antibody [RA3.3A 1/6.1] {PE/Cyanine5} (STJA0000664)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	19264 5788
Gene Symbol	Ptprc PTPRC
UniProt ID	PTPRC_MOUSE PTPRC_HUMAN

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine splenocytes were stained with Anti-CD45R/B220 antibody (STJA0000664) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.