

Anti-CD183/Cxcr3 antibody [CXCR3-173] {PE/Cyanine5} (STJA0000679)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	12766
Gene Symbol	Cxcr3
UniProt ID	CXCR3_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine splenocytes were stained with Anti-CD183/CXCR3 antibody (STJA0000679) and Anti-Mouse CD3 Monoclonal Antibody (AF488 Conjugated) (left). Splenocytes stained with Anti-Mouse CD3 Monoclonal Antibody (AF488 Conjugated) (right) were used as control.