

Anti-CD49b/Igta2 antibody [DX5] {APC} (STJA0000699)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	APC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	16398
Gene Symbol	Itga2
UniProt ID	ITA2_MOUSE

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Mouse bone marrow cells were stained with Anti-CD49b antibody (STJA0000699) [Used at 0.02 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained bone marrow cells (blank black histogram) were used as control.