

Anti-Ly-6G/Ly-6C antibody [RB6-8C5] {PerCP/Cyanine5.5} (STJA0000724)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	17067 546644 100041546
Gene Symbol	Ly6c1 Ly6g LY6C1_MOUSE LY6G_MOUSE LY6C2_MOUSE <

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine bone marrow cells were stained with Anti-Mouse Gr-1 Monoclonal Antibody (PerCP/Cyanine5.5 Conjugated) (filled gray histogram). Unstained bone marrow cells (empty black histogram) were used as control.