

Anti-Ly-6G/Ly-6C antibody [RB6-8C5] {PE} (STJA0000728)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	17067 100041546 546644
Gene Symbol	Ly6c1 Ly6c2 LY6C1_MOUSE LY6C2_MOUSE LY6G_MOUSE <

ADDITIONAL INFORMATION

Note	FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.
------	---



C57BL/6 murine bone marrow cells were stained with Anti-Mouse Gr-1 Monoclonal Antibody (PE Conjugated)[Used at 0.02 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained bone marrow cells (empty black histogram) were used as control.