

Anti-Ly6C antibody [Mons 1] {FITC} (STJA0000746)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	17067 100041546
Gene Symbol	Ly6c1 Ly6c2
UniProt ID	LY6C1_MOUSE LY6C2_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine bone marrow cells were stained with Anti-Ly6C antibody (STJA0000746) [Used at 0.2 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained bone marrow cells (empty black histogram) were used as control.