

Anti-Ly6C antibody [Monts 1] {PercP} (STJA0000749)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PercP
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	100041546 17067
Gene Symbol	Ly6c2 Ly6c1
UniProt ID	LY6C2_MOUSE LY6C1_MOUSE

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



C57BL/6 murine bone marrow cells were stained with Anti-Ly6C antibody (STJA0000749) [Used at 0.2 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained bone marrow cells (empty black histogram) were used as control.