

## Anti-CD45.2/Ptprc antibody [104.2] {FITC} (STJA0000764)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Mouse

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">19264</a>
Gene Symbol	<a href="#">Ptprc</a>
UniProt ID	<a href="#">PTPRC_MOUSE</a>

### ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



C57BL/6 murine splenocytes were stained with Anti-CD45.2 antibody (STJA0000764) [Used at 0.2  $\mu$ g/ $10^6$  cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.