

Anti-CD1A antibody [OKT-6] {PerCP/Cyanine5.5} (STJA0000803)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	909
Gene Symbol	CD1A
UniProt ID	CD1A_HUMAN

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Human T lymphoblastic leukemia cell line Molt-4 were stained with Anti-CD1a antibody (STJA0000803) (filled gray histogram) or Mouse IgG1 Isotype Control PerCP/Cy5.5 (empty black histogram).