

Anti-CD279/Pdcd1 antibody [29F.1A12] {PE/Cyanine5.5} (STJA0000824)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES


| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 5 Mu L/Test |
| Conjugation | PE/Cyanine5.5 |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 18566 |
| Gene Symbol | Pdcd1 |
| UniProt ID | PDCD1_MOUSE |

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Con-A stimulated C57BL/6 splenocytes (3 days) were stained with Anti-CD279/PD-1 antibody (STJA0000824) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.