

## Anti-CD279/Pdcd1 antibody [29F.1A12] {AF647} (STJA0000826)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	AF647
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">18566</a>
Gene Symbol	<a href="#">Pdcd1</a>
UniProt ID	<a href="#">PDCD1_MOUSE</a>

### ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Con-A stimulated C57BL/6 splenocytes (3 days) were stained with Anti-CD279/PD-1 antibody (STJA0000826) (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.