

Anti-CD279/Pdcd1 antibody [29F.1A12] {PE/Cyanine5.5} (STJA0000831)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	PE/Cyanine5.5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	18566
Gene Symbol	Pdcd1
UniProt ID	PDCD1_MOUSE

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Con-A stimulated C57BL/6 splenocytes (3 days) were stained with Anti-CD279/PD-1 antibody (STJA0000831) [Used at 0.2 μ g/ 10^6 cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.