

## Anti-Cd274/PD-L1 antibody [10F.9G2] {APC} (STJA0000838)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	0.2 mg/mL
Conjugation	APC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">60533</a>
Gene Symbol	<a href="#">Cd274</a>
UniProt ID	<a href="#">PD1L1_MOUSE</a>

### ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



C57BL/6 murine splenocytes were stained with Anti-Mouse PD-L1 Monoclonal Antibody (APC Conjugated)[Used at 0.05  $\mu$ g/ $10^6$  cells dilution] (filled gray histogram). Unstained splenocytes (empty black histogram) were used as control.