

Anti-Cd274/PD-L1 antibody [29E.2A3] {PE} (STJA0000840)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	29126
Gene Symbol	CD274
UniProt ID	PD1L1_HUMAN

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



PHA-stimulated (3 days) human peripheral blood lymphocytes were stained with Anti-CD274/PD-L1 antibody (STJA0000840) (filled gray histogram) or Mouse IgG1 Isotype Control PE (empty black histogram).