

Anti-CD206/Mrc1 antibody [C068C2] {APC} (STJA0000843)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ICFC |
| Host / Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES


| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 5 Mu L/Test |
| Conjugation | APC |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 17533 |
| Gene Symbol | Mrc1 |
| UniProt ID | MRC1_MOUSE |

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Mouse abdominal macrophages elicited by starch broth were stained with Anti-CD206 antibody (STJA0000843) (filled gray histogram). Unstained macrophages (blank black histogram) were used as control.