

Anti-CD206/Mrc1 antibody [C068C2] {APC} (STJA0000844)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ICFC |
| Host / Source | Rat |
| Reactivity | Mouse |

PRODUCT PROPERTIES


| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | 0.2 mg/mL |
| Conjugation | APC |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Gene ID | 17533 |
| Gene Symbol | Mrc1 |
| UniProt ID | MRC1_MOUSE |

ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Mouse abdominal macrophages elicited by starch broth were stained with Anti-CD206 antibody (STJA0000844) (filled gray histogram). Unstained macrophages (blank black histogram) were used as control.