

Anti-CD69 antibody [FN50] {PE/Cyanine7} (STJA0000877)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine7
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	969
Gene Symbol	CD69
UniProt ID	CD69_HUMAN

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



PMA and ionomycin-stimulated (4h) Jurkat cells were stained with Anti-CD69 antibody (STJA0000877) (filled gray histogram) or Mouse IgG1, κ Isotype Control (PE/Cyanine7 Conjugated) (empty black histogram).