

## Anti-CD45RO/PTPRC antibody [UCHL1] {FITC} (STJA0000881)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">5788</a>
Gene Symbol	<a href="#">PTPRC</a>
UniProt ID	<a href="#">PTPRC_HUMAN</a>

### ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Human peripheral blood lymphocytes are stained with Anti-CD45RO antibody (STJA0000881) (filled gray histogram) or Mouse IgG2a, κ Isotype Control (FITC Conjugated) (empty black histogram).