

Anti-CD122/IL-3RB antibody [TU27] {FITC} (STJA0000933)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	FITC
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	3560
Gene Symbol	IL2RB
UniProt ID	IL2RB_HUMAN

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Human peripheral blood lymphocytes were stained with Anti-Human CD56 Monoclonal Antibody (APC Conjugated) and Anti-Human CD122 Monoclonal Antibody (FITC Conjugated).