

## Anti-CD197/CCR7 antibody [G043H7] {PE/Cyanine7} (STJA0000946)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine7
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">1236</a>
Gene Symbol	<a href="#">CCR7</a>
UniProt ID	<a href="#">CCR7_HUMAN</a>

### ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Human peripheral blood lymphocytes are restained with Anti-Human CD197 Monoclonal Antibody (PE/Cyanine7 Conjugated) and Anti-Human CD3 Monoclonal Antibody (FITC Conjugated)