

Anti-CD90/THY1 antibody [5E10] {PerCP/Cyanine5.5} (STJA0000961)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PerCP/Cyanine5.5
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

TARGET INFORMATION

Gene ID	7070
Gene Symbol	THY1
UniProt ID	THY1_HUMAN

ADDITIONAL INFORMATION

Note FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.



Human erythroleukemic cell line (HEL) were stained with Anti-CD90 antibody (STJA0000961) (filled gray histogram) or Mouse IgG1, κ Isotype Control (PerCP/Cyanine5.5 Conjugated) (empty black histogram).