

## Anti-Cd14 antibody [Sa14-2] {PE} (STJA0000979)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Rat
Reactivity	Mouse

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	5 Mu L/Test
Conjugation	PE
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	Store between 2-8°C and protecting from prolonged exposure to light. Do not freeze this product.

### TARGET INFORMATION

Gene ID	<a href="#">12475</a>
Gene Symbol	<a href="#">Cd14</a>
UniProt ID	<a href="#">CD14_MOUSE</a>

### ADDITIONAL INFORMATION

Note **FOR SCIENTIFIC EDUCATIONAL RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR OTHER MEDICAL APPLICATIONS.**



Mouse abdominal macrophages elicited by starch broth were stained with Anti-CD14 antibody (STJA0000979) (filled gray histogram). Unstained macrophages (blank black histogram) were used as control.